



OPTI-farrow II

OPTI-farrow II is a farrowing pen developed for loose sows and with optimal protection of the newborn pigs. The pen is designed according to the latest German requirements. The farrowing pen is equipped with a farrowing crate to control the sow's movements – typically the first three to four days around farrowing. The crate can be opened for greater freedom of movement for the sow. It is placed diagonally in the pen and is designed for the space requirements of today's large sows while providing amply space for nursing large litters. The liftable front plate in the creep area makes it easy to treat and protect piglets from the sow's movements. The creep area is positioned so that the cover and front plate can be operated from the walkway. This location makes it easier to monitor and treat piglets without interference from the sow. In addition, biosecurity is improved when the staff does not need to enter the pen.

**We call ourselves
pig experts and
we stand by that**



Open crate



Farrowing section



Closed crate



Open crate

Facts

- Pen length: 240-330 cm
- Pen width: 240-280 cm
- Height of penning: 50/100 cm
- Length of crate from trough: 220 cm
- Tilttable bars with noise reduction
- Creep area with two-part cover

- Sow trough V2A 21 L
- Water system in corrosion resistant materials
- High yielding sow water nipple minimum 10 l/min
- Water system for piglets integrated in penning
- Water nipple or water/milk cup for piglets

Vissing Agro A/S

Energivej 5, DK-8740 Brødstrup
www.vissingagro.dk



download
our
brochures

